

LOWER



CORE



UPPER



ARM CURL

Designed for bicep isolation, the Selectorized Ease Line Arm Curl provides a stable and effective workout with contoured pads and adjustable resistance.

FREE—FLOATING HANDLES

Comfortably accommodate to different forearm lengths.

ANATOMICAL ARM PAD


Isolated & corrects muscle imbalances and adds training variety.

GAS—ASSISTED SEAT

Adjustment with 4—Bar linkage design provides easy and secure adjustability.

VARIABLE CAM RESISTANCE


Mimics natural strength of the biceps for better muscle activation and stimulation.



Product Weight	217kg
Footprint (L x B x H)	139 x 95 x 150cm
Weight Stack	109kg

CHEST PRESS

Prioritize upper body strength with the Selectorized Ease Line Chest Press, featuring biomechanically sound design for effective chest engagement and variable weight options.



CONVERGING ARC OF MOVEMENT

Allows for greater overall muscle development. Greater peak contraction and top and max stretch at the bottom.

INDEPENDENT ARMS

Isolates & corrects muscle imbalances and adds training variety.

ADJUSTABLE BACK PAD

Can be moved forward or back to achieve full range of motion.

Product Weight	242kg
Footprint (L x B x H)	111 x 143 x 150cm
Weight Stack	109kg



LAT PULLDOWN


Our Selectorized Ease Line Lat Pulldown station offers a compound upper body workout, with adjustable resistance and ergonomic seating for enhanced muscle development.

INDEPENDENT ARM

Isolated & corrects muscle imbalances and adds training variety.

ANCHOR THIGH PAD


Keeps exercisers fixed and stable throughout the movement.



Product Weight	260kg
Footprint (L x B x H)	152 x 143 x 185cm
Weight Stack	109kg

SHOULDER PRESS

The Selectorized Ease Line Shoulder Press features a user—friendly setup to target deltoids effectively, offering a stable workout with customizable resistance levels.



INDEPENDENT ARMS

Isolates & corrects muscle imbalances and adds training variety.

THREE GRIP POSITIONS

Accommodates shoulder width variation and different areas of the deltoid to be targeted.

RECLINED SEAT

A slight lean in the starting position enforces muscle activation and increased stability.

COUNTER BALANCED ARMS


Reduces starting resistance for easy lift—off.



Product Weight	254kg
Footprint (L x B x H)	157 x 153 x 150cm
Weight Stack	109kg

BACK EXTENSION

The Selectorized Ease Line Back Extension is crafted for comfort and effectiveness, focusing on lower back muscles with precision and featuring adjustable resistance for progressive training sessions.



ADJUSTABLE START

Allows a selection for range of motion. ROM adjustment is easily selected with the hand release and accommodates individual preferences or physical limitations.

DUAL FOOT PLATFORMS

Accommodate a wider range of users. Dual foot supports provide support and stabilization during exercise without the need for adjustments.

ERGONOMIC SEAT

Ensures correct positioning on the equipment without the need for adjustments.

BACK PAD

The back pad design encourages extension of the lower back muscles while minimizing hip extension.

Product Weight	243kg
Footprint (L x B x H)	106 x 129 x 150cm
Weight Stack	109kg



ABDOMINAL

Designed for core engagement, the Selectorized Ease Line Abdominal station offers a focused workout to strengthen and tone the midsection with intuitive resistance adjustments.

ERGONOMIC SEAT

Ensures correct positioning on the equipment without the need for adjustments.

DUAL SWIVEL PULLEYS


Support both “rotating” and traditional “in—line” abdominal training.

ERGONOMIC KNEELING POSITION

Stabilizes the user’s body and secures user in place or The dual foot supports accommodate users of every height and prevent the hip flexors from being activated.

BACK PAD


Stabilizes the pelvis for better isolation of the abdominal muscles.



Product Weight	218kg
Footprint (L x B x H)	107 x 118 x 157cm
Weight Stack	109kg

LEG PRESS

Experience the Selectorized Ease Line Leg Press, combining ergonomic excellence with variable resistance, enabling users to enhance leg strength and muscular endurance in a supported position.



ADJUSTABLE SEAT CARRIAGE

Allows starting position adjustment for a wider range of users.

ROBUST LINEAR GUIDE ROD

Ensures smooth, solid movement with low friction.

LOW—LYING DESIGN

For easy entry and exit.

LARGE FOOT PLATFORM

Increase exercise variety.

ANGLE OF RAIL INCLINATION

Keeps knees in a biomechanical position for safer and higher—quality quad activation.




Product Weight	347.2kg
Footprint (L x B x H)	190 x 111 x 150cm
Weight Stack	170kg

Selectorized Ease Line EL—08

Intenza Fitness_____ Strength Equipment

SEATED LEG CURL

The Selectorized Ease Line Seated Leg Curl allows for concentrated hamstring training through an adjustable seat and backrest, enhancing user performance.



SELECTION OF MOTION ADJUSTMENT

Accommodates wider range of users.

ADJUSTABLE THIGH PAD

Assisted for stabilization during exercise.

SELF—PIVOTING TIBIA PAD


Allows the ankle to comfortably move around.



Product Weight	261kg
Footprint (L x B x H)	116 x 166 x 150cm
Weight Stack	109kg

LEG EXTENSION

The Selectorized Ease Line Leg Extension machine delivers targeted quadriceps conditioning with its smooth, guided motion for consistent tension and muscle activation.



SELECTION OF MOTION ADJUSTMENT

Accommodates a wider range of users.

SELF—PIVOTING TIBIA PAD

Makes user operation easy and simple.

RECLINED SEAT AND BACK PAD

Help relax hamstrings for full quadriceps contraction.

VARIABLE CAM


Mimics natural strength of the hamstring for better muscle activation and stimulation.

Product Weight	241kg
Footprint (L x B x H)	116 x 141 x 150cm
Weight Stack	109kg



SEATED ROW

The Selectorized Ease Line Seated Row targets the back, ensuring proper posture and peak engagement, complete with selectable weights for tailored workouts.



VERSATILE GRIP SYSTEM

Ensures targeting lats and rear deltoids. Rotating/pivot arm feature.

INDEPENDENT ARMS

Isolated & corrects muscle imbalances and adds training variety.

ADJUSTABLE CHEST PAD

Accommodates for wider range of users.

ROTATING HANDLES

Enable natural movement, allowing users to perform a comprehensive range of motions that target various muscle groups

FEET PAD


Provides a grippy surface for better stability.

Product Weight	248kg
Footprint (L x B x H)	151 x 120 x 150cm
Weight Stack	109kg



SEATED DIP

Focused on triceps and chest, the Selectorized Ease Line Seated Dip station provides a secure position with adjustable resistance to suit all fitness levels.



ROTATING HANDLES

Accommodates to differing shoulder widths.

4—BAR LINKAGE MOVEMENT STRUCTURE

Deliveries a natural movement with minimized stress on wrist joint.

ANGLED BACK PAD

Provides stability to the user eliminating the need for additional adjustment.



Product Weight	242kg
Footprint (L x B x H)	123 x 107 x 150cm
Weight Stack	109kg

WEIGHT ASSISTED CHIN UP/DIP

The Selectorized Ease Line Weight Assisted Chin Up/Dip station aids in perfecting chin—ups and dips by providing adjustable weight assistance for all user strengths.

MULTIPLE GRIP POSITION

Multiple grip positions on Chin—up increase variety of exercises. Dual handle (neutral and prone) for different muscle activation.

HANDLE DESIGN

Dip handles swivel in narrow and wide positions to target different muscles of training.


STANDING DESIGN

Help align whole body better throughout range of motion when compared to kneeling designs.

STEPS DESIGN

Assist with easy entry and exit.

Product Weight	262kg
Footprint (L x B x H)	196 x 99 x 249cm
Weight Stack	109kg




Selectorized Ease Line EL—13

Intenza Fitness_____ Strength Equipment

HIP ABDUCTION/ADDUCTION

Our Selectorized Ease Line Hip Abduction/Adduction station caters to both inner and outer thigh muscles, providing an effective workout with easy—to—modify resistance settings.



DUAL EXERCISES MODE

For better usage and floor space efficiency.

LOW—LYING DESIGN

For easy entry and exit.

ADJUSTABLE STARTING POSITION

Accommodates a wide range of exercisers.



Product Weight	247kg
Footprint (L x B x H)	159 x 82 x 150cm
Weight Stack	109kg

PEC FLY/REAR DELT

Strengthen and tone with the Selectorized Ease Line Pec Fly/Rear Delt, allowing full range of motion and resistance customization for comprehensive chest and shoulder workouts.

SIMPLE ADJUSTMENT

Simple adjustment to switch between exercises.

ARTICULATING ARMS


Accommodates wide range of users.

DUAL HANDLE (NEUTRAL & PRONE)

For different muscle activation.

VARIABLE CAM RESISTANCE


Ensures a direct and secure link to the weight stack.



Product Weight	232kg
Footprint (L x B x H)	141 x 130 x 207cm
Weight Stack	109kg

PRONE LEG CURL

Featuring an auto—adjusting roller pad, the Prone Leg Curl conforms to users of various sizes, ensuring a comfortable and personalized workout. The machine allows users to maximize hamstring engagement while minimizing stress on the lower back.



LOWER—BACK RELIEF

The divergent hip and chest pad design minimizes lower—back stress.

AUTO—ADJUST PAD

The free—floating roller pad automatically adjusts to accommodate users of different sizes without the need for manual adjustment.

NATURAL STRENGTH CAM

The variable cam is designed to replicate the user's natural strength curve throughout the motion.

SECURE GRIP HANDLES

Angled rubber grip handles with aluminum fittings are strategically designed to secure the user in place during exercise.

SINGLE ADJUSTMENT

Designed for ease of use with a single adjustment feature.

Product Weight	219kg
Footprint (L x B x H)	156 x 111 x 150cm
Weight Stack	109kg



TRICEPS EXTENSION

Designed to provide a comfortable and biomechanically sound workout for the triceps, the Ease Line Tricep Extension offers a superior combination of ergonomic design, adjustability, and durability. Its biomechanically sound and ergonomic design effectively isolates and activates the triceps.

COMFORT—FIRST ERGONOMICS


Minimize strain on the joints and reduce the risk of injury with the ergonomic press pad.

EASY ADJUSTMENTS

Allow users of various heights and sizes to find their optimal positioning with easily adjustable seat and arm pads.

BIOMECHANICALLY SOUND

Ensure proper form throughout the range of motion, promoting effective tricep isolation and better muscle activation.



Product Weight	204.1 kg
Footprint (L x B x H)	114 x 95 x 150cm
Weight Stack	109kg